









ECG ELECTRODES - ADULT FOAM			
Gel Type	Adhesive Type	Metal Snap	Radiotranslucent (RTL) Carbon Snap
Conductive Adhesive Gel	General Purpose	 <p>1800 Suretrace &amp; 3500 Profile Foam</p>	 <p>1800C Suretrace RTL &amp; 2100-100 E-Trode</p>
	Diaphoretic/Diagnostic	 <p>1870 Positrace</p>	 <p>1870C Positrace RTL</p>
Wet Gel	General Purpose	 <p>1510 Exactrace</p>	 <p>1500 Dyna/Trace</p>
	Diaphoretic/Diagnostic	 <p>1590 Instatrace</p>	 <p>1550 Dyna/Trace Stress &amp; 1555 Dyna/Trace Mini</p>