

DERM 101

Bovie[®]



MAKING THE OFFICE PRACTITIONER MORE PRODUCTIVE



HIGH-FREQUENCY DESICCATOR

ECONOMICAL

Affordably priced to minimize the facility's expense and maximize their ROI.

AMPLE POWER

Ten (10) watts of monopolar power which are adjustable in 1/10th watt increments.

SIMPLICITY

The thoughtfully designed unit allows Practitioners the ability to perform minor skin procedures with an easy-to-operate device, to increase their practice revenue.

EQUIPPED

Packaged complete and ready to use. By the time the practitioners utilize all 18 electrodes they will have paid for the unit, twice.

RELIABILITY

The Derm 101 is backed by Bovie's standard four (4) year warranty.


DERM 101
INCLUDES

			
DERM 101 Unit	A902 3-Button Handpiece	A804 Sharp Dermal Tips, Non-Sterile (5ea)	A805 Sharp Dermal Tips, Sterile (2ea)
			
A806 Blunt Dermal Tips, Non-Sterile (5ea)	A807 Blunt Dermal Tips, Sterile (2ea)	A806DE Derm-Elite™ Premium Blunt Dermal Tip, Non-Sterile (2ea)	A807DE Derm-Elite™ Premium Blunt Dermal Tip, Sterile (2ea)
			
A910 Disposable Handpiece Sheath, Non Sterile (2ea)	A837 Wall Mount Kit (1ea)	Medical Grade Power Cord	User Manual (CD Format)

DERM 102
& ACCESSORIES

The Bovie® Derm 102 has an additional Bipolar mode offering up to 10 watts of power in 1/10th watt increments.





		
A803 Footswitch*	Reusable Bipolar Forceps* †	A827V Bipolar Forcep Cable*

*Sold Separately † More styles available

Additional accessories available—contact your Bovie representative for details

DERM 101 and DERM 102 Technical Specifications...

Dimensions and Weight:	Mains Connection:
Height 2.7" (69mm) Width 7.2" (184mm) Length 8.2" (208mm) Weight < 1.4 kg (< 3 lbs)	Main Voltage: 100 – 240 VAC ± 10% Main Frequency: 50 – 60 Hertz Main Current: 320mA Max. Power Consumption: 28 VA
Safety:	 
Basic Construction: In accordance with EN 60601-1 Mode of Operation: Intermittent operation Output Type: TYPE BF Duty Cycle: 10sec on / 30sec off Main Fuses: T 400mAL, 250V	

TECHNICAL
SPECIFICATIONS



Bovie Medical Corporation • 5115 Ulmerton Road • Clearwater, FL 33760 • www.boviemed.com • sales@boviemed.com
U.S. Phone 1-800-537-2790 • Fax 1-800-323-1640 • Int'l Phone +1-727-384-2323 • Fax +1-727-347-9144

Bovie® is a registered trademark of Bovie Medical Corporation